

UpDate!

1896 SCHWERT
SWORD SYNONYM FÜR QUALITÄT & VERTRAUEN
SYNONYM FOR QUALITY & TRUST

125
YEARS
1896-2021

SCHWERT Komposite Fissurenzieher im **COLOR** CODE Griff SCHWERT composite cone burnisher in **COLOR** CODE handle

Modern Therapy · Modern Instruments



SCHWERT
SWORD · ESPADA
EPEE · SPADA

MADE IN GERMANY

schwert.com

i

ULTRA Komposite · ULTRA composite

Die schwarze ULTRA Beschichtung ist im Vergleich zu Standardinstrumenten extrem hart und kratzbeständig. Die blendfreie und glatte Oberflächenstruktur verhindert störende Lichtreflexionen und ein Anhaften von Kompositen. Die chemischen Eigenschaften der Beschichtung führen bei der Anwendung von Lösungsmitteln oder Sterilisation zu keinen unerwünschten Reaktionen.

The black ULTRA coating is extremely hard and scratch resistant compared to standard instruments. The smooth and anti glare surface avoids disturbing reflections of light and is less sticky for composite materials. The coating does not react with cleaning solutions or throughout the sterilization process, due to its chemical characteristics.

ULTRA



4085-12Ti

Komposite Fissurenzieher
Composite cone burnisher


Color Code Griff grün, auswechselbare Spitzen, REF 4084-121Ti + 4084-131Ti

Ø 1,8 / 2,2 mm, schwarz beschichtet, 17 cm

Color Code green handle, tips exchangeable, REF 4084-121Ti + 4084-131Ti,

Ø 1,8 / 2,2 mm, black coated, 17 cm

Ergänzend · In addition



4085-00Ti

ULTRA Komposite Set, beinhaltet: 4085-01Ti,
4085-02Ti, 4085-03Ti, 4085-04Ti, Waschtray 2599-05

ULTRA Composite set, consisting of: 4085-01Ti,
4085-02Ti, 4085-03Ti, 4085-04Ti, tray 2599-05



3300-00

Instrumentenschlüssel, rostfrei, für auswechselbare Einsätze, ohne Handschutz, 12,5 cm
Instrument wrench, stainless, for exchangeable inserts, without hand guard, 12.5 cm