

UpDate!


1896 SCHWERT
SWORD SYNONYM FÜR QUALITÄT & VERTRAUEN
SYNONYM FOR QUALITY & TRUST

LUCAS Küretten, SCHWERT Design LUCAS Currettes, SCHWERT Design

Modern Therapy · Modern Instruments



SCHWERT
SWORD · ESPADA
EPEE · SPADA



MADE IN GERMANY

schwert.com



SCHWERT
SWORD · ESPADA
EPEE · SPADA

UpDate!

i

Lucas Küretten · LUCAS Currettes

Scharfe Löffel nach Lucas, auch Lucas Schaber oder Lucas Kürette genannt, werden im Regelfall benutzt um Zysten zu entfernen oder zum Ausschaben von z. B. der Alveole nach einer Extraktion.

Die neuen SCHWERT Lucas Küretten sind mit dem bewährten SCHWERT Design Griff ausgestattet, um dem Behandler ein optimales Handling während des Arbeitens zu bieten. Des Weiteren sind die Arbeitenden grazil ausgearbeitet um einen perfekten Zugang zum Operationsgebiet zu erhalten. Außerdem lässt sich die glatte Oberfläche nach der Behandlung optimal reinigen.


Die Lucas Küretten, welche es in den Größen 1,5 mm, 2,5 mm und 3,0 mm gibt, sind nun auch mit gezahnten Arbeitenden verfügbar. Die feine Zahnung ermöglicht dabei eine ideale Behandlung, im Bereich der Zystenentfernung. Durch die Zahnung ist das Instrument wesentlich aggressiver, was ein Wegrutschen, an der zu behandelnden Stelle verhindert.

Sharp spoons after Lucas or so called Lucas currettes are mostly used to remove cysts or to clean the alveole after an extraction.

The new SCHWERT Lucas currettes are made with the proven SCHWERT design handle, to give a perfect handling for the doctor during the treatment. The working tips are made delicate to allow a perfect approach in the area of the treatment. Additionally the very smooth surface is very easy to clean.

These Lucas currettes, which are available in the sizes: 1.5 mm, 2.5 mm and 3.0 mm are now also available with serrated working tips. The fine serration allow a very good treatment, especially in the area of removing cysts. Due to the serration the instrument is a lot more aggressive, which avoids a slipping away from the spot which has to be treated.

ungezahlt · non-serrated



gezahnt · serrated






Abbildung
Illustration




Artikel-Nr.
Item-No.

Artikelbezeichnung
Article Description

Ungezahnte LUCAS Küretten, SCHWERT Design · Non-serrated LUCAS Currettes, SCHWERT Design

	1372-0	Ø 1,5 mm, ungezahlt Ø 1.5 mm, non-serrated
	1372-1	Ø 2,5 mm, ungezahlt Ø 2.5 mm, non-serrated
	1372-2	Ø 3,0 mm, ungezahlt Ø 3.0 mm, non-serrated

Gezahnte LUCAS Küretten, SCHWERT Design · Serrated LUCAS Currettes, SCHWERT Design

	1372-10	Ø 1,5 mm, gezahnt Ø 1.5 mm, serrated
	1372-11	Ø 2,5 mm, gezahnt Ø 2.5 mm, serrated
	1372-12	Ø 3,0 mm, gezahnt Ø 3.0 mm, serrated

